

## SCOUT and CODE

### Download and installation remarks

Aachen, January 2004


Dear customer,

you can download the latest versions of our SCOUT and CODE programs from our homepage [www.mtheiss.com](http://www.mtheiss.com). Go to the 'Support' section and from there to 'Download'. Each program has its own section with a short description of the items you can download.


There are two types of program downloads:

- 1.The full setup comes in a large zip-file. It contains a setup routine that will install the current version of the software like it is distributed on our setup CDs. You will have to de-install any previous version before you execute the new setup. The full setup is not updated as frequently as the small download file described below.
- 2.The small download file is a zip-file as well. It contains the latest program executable and sometimes a few additional files which belong to the latest version of the software. No de-installation is necessary for this download version: Just copy the files to the right destination folders, overwriting existing files of the same name.

Starting with object generation 2.20 we support some Ocean-Optics spectrometers in SCOUT and CODE directly. The built-in Ocean-Optics driver searches for a file called OOIDRV.INI in the Windows directory (or WINNT directory in the case of Windows 2000). If the file is not present the following dialog will show up when you start SCOUT (or CODE) for the first time:



If you don't use Ocean-Optics hardware set the A/D Converter Type to 'NONE' as in the following example



and press OK. SCOUT will save the file OOIDRV.INI and you will never see this dialog again.

With best regards

W. Theiss